

**Please read carefully and keep this package insert for future reference.  
It is recommended that all eye care practitioners  
review this information with their patients.**

### **Biotrue® Hydration Boost Contact Lens Rehydrating Drops**

We believe that the best inspiration comes from nature. That's why Biotrue Contact Lens Drops are inspired by the biology of the eye. Our innovative drops are Preservative Free and are balanced to match the pH of healthy tears. Our formula is infused with naturally inspired ingredients, including hyaluronan (HA), a natural moisturizer found in tears<sup>†</sup>, an antioxidant to protect the HA and an electrolyte found naturally in tears. Biotrue Contact Lens Drops help prevent the formation of deposits and provide instant moisture that lasts up to 8 hours\*. They work well with soft contact lenses, including daily disposables, and rigid gas permeable (RGP) lenses.

#### **DIRECTIONS:**

To lubricate and rewet contact lenses follow these steps:

STEP 1: During lens wear, remove the cap, place 1 or 2 drops in the eye and blink 2-3 times as needed.

STEP 2: If the lens still does not feel comfortable, add another drop.

To help prevent deposit buildup: Place 2 drops in each eye and blink 2-3 times, 4 times a day.

**CONTENTS:** Biotrue Contact Lens Drops is a sterile, preservative free, buffered, isotonic, aqueous solution that contains boric acid, sodium borate, potassium chloride, erythritol, poloxamine, glycerol and hyaluronan.

**ACTIONS:** Biotrue Contact Lens Drops lubricate and rewet soft (hydrophilic) contact lenses including silicone hydrogel contact lenses, and RGP (silicone acrylate [SA] and fluorosilicone acrylate [FSA]) contact lenses. When used 4 times a day, Biotrue Contact Lens Drops helps prevent deposit build-up on lenses.

#### **INDICATIONS (Uses):**

Biotrue Contact Lens Drops are indicated for use to lubricate and rewet soft (hydrophilic) contact lenses, including silicone hydrogel as well as SA and FSA RGP contact lenses during wear. The product may be used with daily or extended wear, and disposable lenses.

#### **CONTRAINDICATIONS (Reasons not to use):**

Do not use if you are allergic to any ingredient in this product.

#### **PROBLEMS WITH CONTACT LENSES AND LENS CARE PRODUCTS COULD RESULT IN CORNEAL INFECTION AND/OR ULCER AND LEAD TO LOSS OF VISION.**

It is essential that you follow your eye care practitioner's directions and all labeling instructions for proper use and care of lenses and lens care products, including the lens case. Eye problems, including corneal ulcers, can develop rapidly and lead to loss of vision. Daily wear lenses are not indicated for overnight wear and should not be worn while sleeping. Clinical studies have shown that the risk of serious adverse reactions is increased when these lenses are worn overnight. Extended wear lenses should be regularly removed for cleaning and disinfection or for disposal and replacement on the schedule prescribed by your eye care practitioner. Clinical studies have shown that there is an increased incidence of serious adverse reactions in extended wear contact lens users as compared to daily wear contact lens users. Studies have also shown that the risk of serious adverse reactions increases the longer extended wear lenses are worn before removal for cleaning and disinfection or for disposal and replacement. Studies have also

shown that smokers have a higher incidence of adverse reactions. If you experience eye discomfort, excessive tearing, vision changes or redness of the eye, immediately remove your lenses and promptly contact your eye care practitioner. It is recommended that you see your eye care practitioner regularly to assure the continued health of your eyes.

#### **DISCARD DATE:**

- Discard any remaining solution within 90 days after opening.

#### **WARNING:**

- Using your contact lens drops beyond the discard date could result in contamination of the solution and can lead to severe infection, vision loss or blindness.

#### **WATER ACTIVITY:**

- Do not expose your contact lenses to water while you are wearing them.
- Ask your eye care practitioner for recommendations about wearing your lenses during any activity involving water.

#### **IMPORTANT SAFETY INFORMATION:**

- Always follow the product directions for use. Failure to follow product directions may lead to vision loss.
- Visit your eye care practitioner regularly.
- Always wash and dry hands before handling lenses.
- Biotrue Contact Lens Drops will not disinfect your lenses.
- To avoid contamination, do not touch tip of container to any surface. Replace cap after using.
- Do not alter the special bottle, the tip opening or the liner inside the cap which maintain sterility without preservatives.

#### **PRECAUTIONS:**

- Replace cap after using.
- Store solution at room temperature.
- Use before the expiration date marked on the carton and bottle.
- Discard any remaining solution within 90 days after opening.
- Keep out of reach of children.

#### **ADVERSE REACTIONS (Problems and What to Do)**

You should be aware that the following problems may occur:

- Eyes stinging, burning, itching (irritation), or other eye pain
- Comfort is less than when lens was first placed on eye
- Abnormal feeling of something in the eye (e.g., foreign body, scratched area)
- Excessive watering (tearing) of the eyes or unusual eye secretions
- Redness of the eyes
- Reduced sharpness of vision (poor visual acuity)
- Blurred vision, rainbows, or halos around objects
- Sensitivity to light (photophobia)
- Dry eyes

#### **IF YOU NOTICE ANY OF THE ABOVE: Immediately remove your lenses.**

- If the discomfort or problem stops, look closely at the lens.
- If the lens is in any way damaged, do not put the lens back on your eye. Place the lens in the lens case and contact your eye care practitioner.
- If the lens has dirt, an eyelash or other foreign body on it, or the problem stops and the lens appears undamaged, thoroughly clean, rinse and disinfect the lens, then reinsert it.

#### **If the problem continues, immediately remove the lens and consult your eye care practitioner.**

If any of the above problems occur, a serious condition such as infection, corneal ulcer, neovascularization or iritis may be present. Seek immediate professional identification of the problem and prompt treatment to avoid serious eye damage.

#### **GOOD LENS CARE HYGIENE:**

- Always wash and dry your hands before handling lenses.
- Clean, rinse and disinfect your lenses each time you remove them.
- Always handle the same lens first, the right or the left, to avoid mix-ups.
- Before first and after each additional use, always empty and rinse lens case and caps with fresh rinsing solution and allow to air dry.
- When wearing disposable lenses, each time you remove a lens, discard it and replace it with a new lens.

**HOW SUPPLIED:** Biotrue Contact Lens Drops are supplied in a 0.33 FL OZ (10 mL) bottle. Bottles and cartons are marked with a lot number and expiration date.

If you are not completely satisfied, need to report a complaint, an adverse event, or need additional product information, call 1-800-553-5340.

Explanation of symbols located at [www.bausch.com/symbols](http://www.bausch.com/symbols)

**Manufactured for:**  
Bausch & Lomb Incorporated  
1400 North Goodman Street  
Rochester, NY 14609 USA  
Biotrue is a trademark of Bausch & Lomb Incorporated or its affiliates.  
© 2023 Bausch & Lomb Incorporated or its affiliates  
[www.biotrue.com](http://www.biotrue.com)  
[www.bausch.com](http://www.bausch.com)

<sup>†</sup>Hyaluronan is sourced from a large scale natural fermentation process.

\*Based on a laboratory study